




SHOWER BAY
by Forward Day



**Shower Safely
Without Remodeling**




SHOWER BAY is a safe and portable shower designed for wheelchair users. Unlike other solutions on the market, Shower Bay offers a true shower experience without requiring dangerous wet-environment transfers or expensive home renovations. The unique design allows for quick assembly in any room of the home.

No Costly Remodel



Shower Bay can be set up in any room of the home.



- ◆ Can be assembled quickly without any tools or need for a contractor or plumber.
- ◆ Water hose extends 10.6m in length for easy access to nearest sink location.
- ◆ Water in-hose easily snaps onto existing sink faucet. Place out-hose in sink or tub basin for drainage. Hoses can be coiled up for storage when not in use.
- ◆ Ramp detaches to create smaller footprint (1140 x 1120mm) between shower uses.



Safe

Ramp allows for safe and smooth access to shower

- ◆ Provides easy ramp access
- ◆ Eliminates dangerous transfers in a wet environment
- ◆ Protects caregiver from physical injury while assisting user



Keep Loved Ones at Home

Enjoy a true shower experience

- ◆ Continue care for a loved one in the comfort of their own home
- ◆ Maintain privacy and dignity behind shower walls
- ◆ Save thousands of dollars on assisted living costs



Shower Bay is a Great Option for:

Home Owners

Those wanting to avoid a costly and time-consuming bathroom remodel

Renters

Those who generally cannot modify or remodel their apartment or home

Aging in Place

Families seeking to avoid the high costs of assisted living facilities and keep their loved one at home

First Floor Living Restrictions

Homes with only a half-bath available on the first floor

Temporary Living Situations

Families that have a loved one with mobility limitations temporarily stay at their home

Injury Rehab Patients

Individuals with the need for a temporary bathing solution during the rehabilitation process

Residential Moves

Users planning to change residences will benefit from the one-time investment in a portable shower

Equipment Rentals

Those looking to rent the shower for time or cost reasons

Testimonials

"It's very thoughtfully designed for user and caregiver needs... I would definitely recommend the Shower Bay to patients and families in need of a non-permanent roll-in shower."

– Allison Joe

Senior Occupational Therapist, Stanford Hospital

"The overall design of Shower Bay is efficient and effective, with easy access and adequate shower space... This product is a very attractive alternative to a full shower remodel, particularly for those with challenging health issues."

– Joe S (ALS patient) and wife, Mary S (caregiver)



Frequently Asked Questions

Where does the water come from?

Shower Bay's unique quick-connect system easily snaps onto most home faucets. The water travels from the faucet through an in-hose to the Shower Bay unit.

Where does the water go?

The built-in pump located at the base of Shower Bay removes the water through the attached out-hose to a nearby sink or tub drain.



Is Shower Bay supposed to fit in my bathroom?

Shower Bay may be used in a bathroom, but it was also designed to work in any room of the home. We understand that most bathrooms are not large enough for adequate wheelchair maneuverability. In most cases Shower Bay is set up in the corner of a bedroom or spare area of the home, such as a den, laundry room or garage.

Do I need a contractor, plumber or permits?

Because there is no installation, a contractor, plumber or building permit is not required.

Does Shower Bay fill up like a bathtub?

Shower Bay is a shower and does not fill up like a bathtub.



View more FAQ's online at www.ShowerBay.com

Features

- ❖ Quick assembly – No tools required
- ❖ Eliminate dangerous wet-environment transfers with roll-in ramp
- ❖ Maintain privacy and preserve dignity behind shower walls
- ❖ Easily snap “quick-connect” adapter onto sink faucet
- ❖ Shower in most rooms of the home using the flexible 35-foot water supply hoses
- ❖ Ramp detaches for compact storage



Specifications

Ramp:	890mm	Length
	1120mm	Width
Shower Base:	1140mm	Length
	1120mm	Width
Shower Base and Ramp:	2032mm	Length
	1120mm	Width
Overall Height:	1270mm	Height
Hose: *	890mm	Length
Power Cord:	300mm	Length

* Longer hose length available upon request

	Shower Bay	Bath Remodel	Bath Bench	Walk-in Tub	Sponge Bath
Safe	✓	✓		✓	✓
Clean	✓	✓	✓	✓	
Easy to Use	✓	✓		✓	✓
Privacy	✓	✓	✓	✓	
No Tools Assembly	✓		✓		✓
Affordable	✓		✓		✓
Immediate Availability	✓		✓		✓
Roll-in Entry	✓	✓			

Dimensions





Australian Distributer:



home modifications

Safety at Home

16-20 Ragless Street, St Marys SA 5042

T: 08 8374 1099

office@hiline.net.au

www.hiline.net.au

Check out more pictures and videos at

www.ShowerBay.com



Made in U.S.A.